INDLAB FUME HOOD

HOTTE INDLAB
INDLAB SORBONNE

INDLAB

The Indlab model was developed to carry out tests with dangerous products, guaranteeing a high reliability in its operation. It always ensures a laminar flow velocity in the inner chamber (working area), avoiding turbulence and the possibility of escaping dangerous gases.

This is the premium model of Industrial Laborum, being an INDLAB Fume Hood certified according to **EN 14175-2** and **EN 14175-3** by TÜV Nord®

Developed according to standards:

EN 14175-1 - Chemistry Fume Hoods - Part 1: Definition and Dimensions

EN 14175-2 - Chemistry Fume Hoods - Part 2: Performance and Safety Requirements

EN 14175-3 - Chemistry Fume Hoods - Part 3: Type Test Methods

EN 13150 - Laboratory Benches - Dimensions, Safety Requirements and Test Methods.

EN 61010-1 - Safety requirements for electrical equipment for measurement,

control and use in the laboratory - Part 1: General requirements

EN 61326-1 - Electrical equipment for measurement, control and use in the laboratory - Electromagnetic compatibility requirements - Part 1: General requirements.

CE marking according to the directives: **2014/30/EU** (Electromagnetic Compatibility) and **2014/35/EU** (low voltage).

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LINE
LINHA
GAMME


This model was developed into a line, which presents different possibilities, according to the normes (EU / US) and user needs, being one of the more complete lines in the market.

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STANDARD
STANDARD
STANDARD


Fume Hood Indlab Standard for a generalized use of chemical material with constant air volume (CAV).

Hotte Indlab Standard para um uso generalizado de material químico, com volume de ar constante (VAC).

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STAINLESS STEEL
AÇO INOX
ACIER INOXYDABLE


Fume Hood Indlab in stainless steel, with constant air volume (VAC), recommended when using radio-chemical substances.

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POLYPROPYLENE
POLIPROPILENO
POLYPROPYLENE


Fume Hood Indlab with interior and worktop in high density Polypropylene, with constant air volume (VAC), advised when using hydrofluoric acids.

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ES.F
ES.F
ES.F


Filtered Fume Hood Indlab with Neutrodine filters, an innovation that allows the absorption of acids, bases and solvents with the same filter.

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U.S. STANDARD
U.S. STANDARD
U.S. STANDARD


Fume Hood Indlab following the norms of the United States of America, with controllers in the towers and services in the interior sides of the Fume Hood.

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ADAPTED
ADAPTADA
ADAPTÉ


Adapted Indlab Fume Hood, with lower surface, controllers in the towers and services on the interior sides of the Fume Hood.

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
Indlab Fume Hood with interior and top in ceramic stoneware, with constant air volume (VAC), advised when there is an intensive use of chemicals.

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
Indlab WALK IN Fume Hood, without work surface and window in the full height of the Fume Hood.

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INDLAB FUME HOOD

High performance chemistry dish designed for use in laboratories that pose high risks to users and the environment.

The INDLAB Fume Hood consists of:

- Exterior body consists of panels made of sheet steel coated with epoxy-polyester resins of high chemical resistance. It serves to support the work surface and the inner body of the hood. In the back space of the bench there is a technical zone where connections to the water, sewage, gas and electricity networks are made;
- Inner housing - Made of various materials with high chemical resistance. It has aerodynamic elements to optimize the flow of air during the extraction.
- Working surface - Made of various materials with high chemical resistance.
- Window - Vertical opening (guillotine type), with a horizontal opening section, allows the user to safely handle dangerous substances inside the Fume Hood.

Our standard models come with a constant air volume controller (CAV), and the services are: a gas spout and two water spouts with the respective controls; also provided with a sink on the top or panels as well as two or three IP55 single-phase sockets. The controls and electrical plugs are placed on the front panels of the underfloor hood, the nozzles being provided inside the hood in panels on the deflectors. There is also the possibility of being supplied with the controls and sockets in the lateral towers, with services on the sides of the interior of the Fume Hood.

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
Air flow
Fluxo de ar
Flux d'air

Standard Fume Hood comes with constant air volume (CAV), however can be equipped with VAV system (variable air volume), which improves the performance of the Fume Hood as well as its overall energy consumption.

In the same way the Fume Hood can also be equipped with automatic window closing system, or with the SCAT chemical waste collection system.

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


Table with 5 columns: Model size, Internal volume, Minimum flow rate, Nominal output diameter, Pressure drop. Rows for models 1200/1300, 1500/1600, 1800/1900.

Standard services
Serviços standard
Services standard


2/3 schuko-type outlets (type F), 230V, protection class IP55**
2/3 tomadas do tipo schuko (tipo F), 230V, classe de proteção IP55**
2/3 sorties de type schuko (type F), 230V, indice de protection IP55 **

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**model 1200 with two sockets and model 1500 and 1800 with 3 sockets. The location of the sockets will be on the front panels of the fume hood, and you must specify if you want in any other location.

One Drip Cup fixed in the service panel or undermounted on the worktop.
Uma pileta fixa ao painel de serviço ou na superfície de trabalho.
Une bénitier fixé dans le panneau de service ou sur le plan de travail.

Electrical information*
Informação eléctrica*
Informations électriques*

The INDLAB Fume Hood is a class I equipment, so all accessible metal parts are connected to the earth terminal, which must be connected to the earth protective circuit of the building.

The installation site must have a circuit-breaker to ensure that in the case of other equipment connected to the circuit-breaker, these shall not be without power in case of maintenance of the Fume Hood.

Table with 5 columns: Model size, Power Supply Voltage, Frequency, Power supply, Installed power. Rows for models 1200/1300, 1500/1600, 1800/1900.

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


STEEL SHEET
CHAPA DE AÇO
TÔLE D'ACIER

COMPACT LAMINATE
COMPACTO FENÓLICO
PHENOLIC COMPACT

COMPACT LAMINATE
COMPACTO FENÓLICO
PHENOLIC COMPACT

STEEL SHEET
CHAPA DE AÇO
TÔLE D'ACIER

COMPACT LAMINATE
COMPACTO FENÓLICO
PHENOLIC COMPACT

POLYPROPYLENE
POLIPROPILENO
POLYPROPYLÈNE

STEEL SHEET
CHAPA DE AÇO
TÔLE D'ACIER

COMPACT LAMINATE
COMPACTO FENÓLICO
PHENOLIC COMPACT

TECHNICAL CERAMIC
GRÉS CERÂMICO
CÉRAMIQUE TECHNIQUE

STEEL SHEET
CHAPA DE AÇO
TÔLE D'ACIER

POLYPROPYLENE
POLIPROPILENO
POLYPROPYLÈNE

TECHNICAL CERAMIC
GRÉS CERÂMICO
CÉRAMIQUE TECHNIQUE

STEEL SHEET
CHAPA DE AÇO
TÔLE D'ACIER

STAINLESS STEEL
AÇO INOX
ACIER INOXYDABLE

STAINLESS STEEL
AÇO INOX
ACIER INOXYDABLE

STEEL SHEET
CHAPA DE AÇO
TÔLE D'ACIER

TECHNICAL CERAMIC
GRÉS CERÂMICO
CÉRAMIQUE TECHNIQUE

TECHNICAL CERAMIC
GRÉS CERÂMICO
CÉRAMIQUE TECHNIQUE

LxDxH*		
1200x1010x2500	105-0007	
1500x1010x2500	105-0008	
1800x1010x2500	105-0009	
1300x1010x2300	105-0028	
1600x1010x2300	105-0029	
1900x1010x2300	105-0030	

1200x1010x2500	105-0001	
1500x1010x2500	105-0002	
1800x1010x2500	105-0003	
1300x1010x2300	105-0022	
1600x1010x2300	105-0023	
1900x1010x2300	105-0024	

1200x1010x2500	105-0004	
1500x1010x2500	105-0005	
1800x1010x2500	105-0006	
1300x1010x2300	105-0025	
1600x1010x2300	105-0026	
1900x1010x2300	105-0027	


1200x1010x2500	105-0013	
1500x1010x2500	105-0014	
1800x1010x2500	105-0015	
1300x1010x2300	105-0034	
1600x1010x2300	105-0035	
1900x1010x2300	105-0036	

1200x1010x2500	105-1001	
1500x1010x2500	105-1002	
1800x1010x2500	105-0049	
1300x1010x2300	105-1021	
1600x1010x2300	105-1022	
1900x1010x2300	105-1023	

1200x1010x2500	105-0043	
1500x1010x2500	105-0044	
1800x1010x2500	105-0045	
1300x1010x2300	105-0064	
1600x1010x2300	105-0065	
1900x1010x2300	105-0066	

L x D x H*		
1200x1010x2500	105-0010	
1500x1010x2500	105-0011	
1800x1010x2500	105-0012	
1300x1010x2300	105-0031	
1600x1010x2300	105-0032	
1900x1010x2300	105-0033	

*Length x Depth x Height
*Comprimento x Profundidade x Altura
*Longueur x Profondeur x Hauteur

STEEL SHEET
CHAPA DE AÇO
TÔLE D'ACIER

FIBERGLASS
FIBRA DE VIDRO
FIBRE DE VERRE

TECHNICAL CERAMIC
GRÉS CERÂMICO
CÉRAMIQUE TECHNIQUE

STEEL SHEET
CHAPA DE AÇO
TÔLE D'ACIER

STEEL SHEET
CHAPA DE AÇO
TÔLE D'ACIER

COMPACT LAMINATE
COMPACTO FENÓLICO
PHENOLIC COMPACT

STEEL SHEET
CHAPA DE AÇO
TÔLE D'ACIER

STEEL SHEET
CHAPA DE AÇO
TÔLE D'ACIER

TECHNICAL CERAMIC
GRÉS CERÂMICO
CÉRAMIQUE TECHNIQUE

STEEL SHEET
CHAPA DE AÇO
TÔLE D'ACIER

CERAMIC TILE
MUSAICO CERÂMICO
MOSAIQUE EN CÉRAMIQUE

TECHNICAL CERAMIC
GRÉS CERÂMICO
CÉRAMIQUE TECHNIQUE

L x D x H*		
1200x1010x2500	105-0001	
1800x1010x2500	105-0003	
1600x1010x2300	105-0023	

MELAMINE
MELAMINA
MÉLAMINE

POLYPROPYLENE
POLIPROPILENO
POLYPROPYLÈNE

TECHNICAL CERAMIC
GRÉS CERÂMICO
CÉRAMIQUE TECHNIQUE

STAINLESS STEEL
AÇO INOX
ACIER INOXYDABLE

STAINLESS STEEL
AÇO INOX
ACIER INOXYDABLE

STAINLESS STEEL
AÇO INOX
ACIER INOXYDABLE


L x D x H*		
1200x1010x2500	105-0016	
1500x1010x2500	105-0037	
1800x1010x2500	105-0017	
1300x1010x2300	105-0038	
1600x1010x2300	105-0018	
1900x1010x2300	105-0039	

1200x1010x2500	105-0052	
1800x1010x2500	105-0053	
1600x1010x2300	105-0054	

1200x1010x2500	105-0055	
1800x1010x2500	105-0056	
1600x1010x2300	105-0057	

1200x1010x2500	105-0046	
1800x1010x2500	105-0047	
1600x1010x2300	105-0048	

L x D x H*		
1200x1010x2500	105-0001	
1800x1010x2500	105-0003	
1600x1010x2300	105-0023	


L x D x H*		
1200x1010x2500	105-0019	
1500x1010x2500	105-0020	
1800x1010x2500	105-0021	
1300x1010x2300	105-0040	
1600x1010x2300	105-0041	
1900x1010x2300	105-0042	

*Length x Depth x Height
*Comprimento x Profundidade x Altura
*Longueur x Profondeur x Hauteur

STANDARD
STANDARD
STANDARD


Standard Indlab Fume Hood for a generalized use of chemical material; exterior in cold-rolled steel sheet, inner lining and worktop in different materials. For a more intensive use of chemical materials it is recommended that the interior and the worktop be in ceramic stoneware, for a better performance of the fume hood over time.

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FLAMMABLES
INFLAMÁVEIS
INFLAMMABLES


Indlab Fume Hood for flammable chemical materials; exterior in cold-rolled steel sheet, inner lining and worktop in stainless steel. It's a Fume Hood advised for those who intend to work with materials from petroleum and some explosive powders, being a Fume Hood that withstands high temperatures while maintaining the safety of its user.

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STRONG ACIDS
ÁCIDOS FORTES
ACIDES FORTS


This Indlab Fume Hood has been developed for acids chemical materials, with exterior in cold-rolled steel sheet, interior and worktop in polypropylene. It's a fume hood advised for those who intend to work with strong acids and very corrosive materials, ensuring stability and safety for the users. However due to the nature of the polypropylene, it is not advisable to use material or equipment with high temperatures, since it will damage the polypropylene.

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RADIOISOTOPE
RADIOISÓTOPO
RADIO-ISOTOPE

This Indlab Fume Hood was developed specifically for radioisotope chemical materials, with exterior, interior and work surface in stainless steel. It's a fume hood advised for those who intend to work with radioisotope, ensuring stability and safety when used. It is also a fume hood advised for those who work in high temperature environments and who need to be in constant cleaning, however, due to the nature of the material, it is not advised to work with some acids.

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