INDLAB
FUME HOOD
HOTTE INDLAB
INDLAB SORBONNE

INDLAB ADAPTED


INDLAB ADAPTADA
INDLAB ADAPTÉE

Model adapted for use by
people with reduced mobility.

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ADAPTED FUME HOOD

A high-performance chemistry Fume Hood suite designed for users of reduced mobility and developed for use in laboratory work that poses a high risk to users and the environment.

The Adapted Fume Hood consists of:

- Exterior body - consists of panels made of sheet steel coated with epoxy-polyester resins of high chemical resistance. It serves to support the work surface and the inner body of the fume hood. In the back space of the bench there is a technical zone where the connections to the water, sewage, gas and electricity networks are made;
- Inner housing - Made of various materials with high chemical resistance. It has aerodynamic elements to optimize the flow of air during the extraction.
- Working surface - Made of various materials with high chemical resistance.
- Sash window - Vertical opening (guillotine type), with a horizontal opening section, allows the user to safely handle dangerous substances inside the fume hood.

Our standard models come with a constant air volume (CAV) controller, and the services are: a gas and two water valves with the respective controls it is also provided with a cup sink in the lateral panels as well as two or three single-phase sockets, IP55. The valves controls and electrical plugs are placed in the uprights of the fume hood, the valves being provided on lateral panels in the interior of fume hood. This fume hood is not suitable for cabinets under the worktop, so the cabinets advised are tall and are next to the fume hood.

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

Air flow


Fluxo de ar
Flux d'air

The Adapted Fume Hood comes standard with constant air volume (CAV), however can be equipped with VAV system (variable air volume), which improves the performance of the Fume Hood as well as its overall energy consumption. In the same way the Fume Hood can also be equipped with automatic window closing system, or with the SCAT chemical waste collection system.

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

Standard services

Serviços standard
Services standard

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2/3 prises de type schuko (type F), 230V, degré de protection IP55 **

*Other services are possible upon request

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Two cold water valves and one natural gas valve, with the respective taps.

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**model 1200 with two sockets and model 1500 and 1800 with 3 sockets. The location of the sockets will be on the front panels of the fume hood, and you must specify if you want in any other location.

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One Drip Cup fixed in the lateral services panel.

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Electrical information

Informação eléctrica
Informations électriques

The Adapted Fume Hood is a class I equipment, so all accessible metal parts are connected to the earth terminal, which must be connected to the earth protective circuit of the building. The terminal strips for electrically connecting the Fume Hood are located in the electrical cabinet on the upright.

The equipment has all necessary magnetothermic protections, namely:

- Thermal protection of the fan(s);
- Protection circuit breaker of the circuit of the sockets;
- Circuit breaker protection of the lighting system;
- Protection circuit breaker for the control circuit.

The installation site must have a circuit-breaker to ensure that in the case of other equipment connected to the circuit-breaker, these shall not be without power in case of maintenance of the Fume Hood.

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

STEEL SHEET
CHAPA DE AÇO
TÔLE D'ACIER

COMPACT LAMINATE
COMPACTO FENÓLICO
PHENOLIC COMPACT

STEEL SHEET
CHAPA DE AÇO
TÔLE D'ACIER

COMPACT LAMINATE
COMPACTO FENÓLICO
PHENOLIC COMPACT

STEEL SHEET
CHAPA DE AÇO
TÔLE D'ACIER

COMPACT LAMINATE
COMPACTO FENÓLICO
PHENOLIC COMPACT

STEEL SHEET
CHAPA DE AÇO
TÔLE D'ACIER

POLYPROPYLENE
POLIPROPILENO
POLYPROPYLENE

STEEL SHEET
CHAPA DE AÇO
TÔLE D'ACIER

STAINLESS STEEL
AÇO INOX
ACIER INOXYDABLE

STEEL SHEET
CHAPA DE AÇO
TÔLE D'ACIER

TECHNICAL CERAMIC
GRÉS CERÂMICO
CÉRAMIQUE TECHNIQUE

STEEL SHEET
CHAPA DE AÇO
TÔLE D'ACIER

POLYPROPYLENE
POLIPROPILENO
POLYPROPYLENE

STAINLESS STEEL
AÇO INOX
ACIER INOXYDABLE

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