ADDITIONAL SYSTEMS

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The automatic sash window system is a window control system, it is a window automation that allows when allied to a VAV extraction system a huge reduction in the cost operation of the fume hood, because, as soon as it is no longer necessary or after some time of the window is fully open with high extraction rates, it automatically lowers and reduces the extraction flow rate.

This system also allows that with just a slight touch up or down the window goes to the maximum or minimum height respectively. All this is automatic thanks to a system that involves various sensors, user presence sensors, objects, as well as height sensors for the window, central control box that compresses all readings and displays the data to the user in the controller; finally we have the engine that accomplishes all these actions. The controller allows you to parameterize the actions that you want the window to perform for different situations.

The automatic window system is something that for a laboratory, that values the safety of its users and for the security of the space, must take into account. It is a system that prevents human errors that can sometimes happen if, for example, leaving the window open while an experience is happening, may contaminate the environment putting users in danger.

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

iberica
Industrial Laborum
lab integrated systems

AIRFLOW SYSTEM

SISTEMA AIRFLOW

SYSTÈME AIRFLOW

AIRFLOW SYSTEM

The integration of the airflow system into the Indlab and Indlab ENS fans allows up to a certain amount of the air to be extracted by the fan to be replaced in the same room by the same fan. This system is a complement to the hoods, but it generates differences in the structure and composition of the hood, so when this system is used, it is considered by Industrial Laborum as a fume hood differentiated from the others, in the case of the INDLAB AIRFLOW FUME HOOD.

This fume hood is certified by TÜV Rheinland, as mentioned above is the incorporation of a system of compensation of the extraction of air, understand whether by compensation air insufflation. This system inflates the air in the room from which the hood is extracting, even, insufflation inside the hood to improve the extraction performance.

In this complement to the hood the volume of inflation is constant, but it is a complement that can be changed to be a system of variable air volume through conventional VAV or through fans with variable speed drives. This system and this hood will remove the HVAC from the lab some effort to monitor the conditions under which the hoods are operated and that can change at any time by putting the HVAC system in some effort.

Airflow Fume Cupboard is certified according to **EN 14175-2:2003** and **EN 14175-3:2003** by TÜV Rheinland®.

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

Scheme of extraction and insufflation

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KITS FOR FUME HOODS

Here we will have all the kits that we can add to our fume hood so that we can work in a much safer environment and able to satisfy our needs. In this chapter we will find the Atex kit, chemical waste kit (S.C.A.T.), gas scrubbers and water neutralizer and finally the fire-trace kit.

The atex kit and fire-trace kit are intended to increase user safety. The atex kit is designed for fume hoods that will work in potentially explosive environments, consisting of several electrical protection equipment. The fire-trace kit consists of fire extinguishing equipment that may occur inside the fume hood. It uses an exponential temperature increase recognition system that triggers a chemical powder extinguisher that will expel on the top of the fume hood to the inside extinguishing the fire.

Gas scrubbers and water neutralizers consist of extra hotte equipment that perform a flushing of the gases, with filtration resulting in a cleaning of the gases of the chemicals present. The water neutralizer balances the pH of the water being poured, whether they are more acidic or alkaline.

The chemical waste system, (S.C.A.T.) is a system that safely allows chemical spills in the surface of the fume hood to a container through a funnel that prevents gases from being expelled and does not imply contamination of water.

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


System that washes the inside of the Fume Hood. Its action allows the automatic washing of the walls and work surface, what constitutes an improvement in the maintenance of the equipment. Not advisable for fume hoods with inner lining with less resistance to water.

Cod. 112-4002

S.C.A.T. system
Sistema S.C.A.T.
Système S.C.A.T.

The S.C.A.T. is a complex system that can be decomposed by the materials in the graphic below. The complete system requires one of each of the components below.


This complete system comes with Signal box that allows to read up two container and the state of capacity in which they are. When one of the containers is full the Signa box gives a loud and luminous warning to the user, thus avoiding the mistake of filling the container too much, where contamination may occur.

All components in contact with chemicals are produced from high chemical resistance polymers.

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


Cod. 031-XXXX

Kit Atex
Kit Atex
Kit Atex

The Atex Kit consists of a set of explosion-proof material for the different fume hoods. This material is for ATEX zone 1, zone 2, zone 21 and zone 22. These components include a waterproof luminaire with two fluoresent lamps and a control button for ventilation and lighting.

For ATEX zones rating 0 or 20 it is recommended that there are no electrical components, with only one control button and the luminaire, passing all other components to the outside of the room. This is because zones 0 and 20 are high-dwelling areas in potentially explosive atmospheres.

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