INDLAB ENS FUME HOOD

HOTTE INDLAB ENS
INDLAB ENS SORBONNE

The Indlab ENS model was developed to carry out tests with dangerous products, guaranteeing a high reliability in its operation. It always ensures a laminar flow inside the work area, avoiding the turbulence and the possibility of escape of dangerous gases.

This Industrial Laborum model is an INDLAB ENS Fume Hood developed according to **EN 14175-2** and **EN 14175-3** standards.

Developed according to standards:

EN 14175-1 - Chemistry Fume Hood - Part 1: Definition and Dimensions

EN 14175-2 - Chemistry Fume Hood - Part 2: Performance and Safety Requirements

EN 14175-3 - Chemistry Fume Hood - Part 3: Type Test Methods

EN 13150 - Laboratory Benches - Dimensions, Safety Requirements and Test Methods.

EN 61100-1 - Safety requirements for electrical equipment for measurement, control and use in the laboratory - Part 1: General requirements

EN 61326-1 - Electrical equipment for measurement, control and use in the laboratory - Electromagnetic compatibility requirements - Part 1: General requirements.

CE marking according to the directives: **2014/30/EU** (Electromagnetic Compatibility) and **2014/35/EU** (low voltage).

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
This model was developed into a line, which presents different possibilities, according to the norms (EU/US) and user needs, being one of the more complete lines in the market.

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STANDARD
STANDARD
STANDARD


Fume Hood Indlab ENS Standard for a generalized use of chemical material with constant air volume (CAV).

Hotte Indlab ENS Standard para um uso generalizado de material químico, com volume de ar constante (VAC).

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STAINLESS STEEL
AÇO INOX
ACIER INOXYDABLE


Fume Hood Indlab ENS in stainless steel, with constant air volume (VAC), recommended when using radio-chemical substances.

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
POLYPROPYLENE
POLIPROPILENO
POLYPROPYLENE

Fume Hood Indlab ENS with interior and worktop in high density Polypropylene, with constant air volume (VAC), advised when using hydrofluoric acids.

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INDLAB ENS FUMEHOOD

Chemistry Fume Hood, designed for use in laboratories that pose high risks to users and the environment.

The Indlab ENS Fume Hood consists of:

- External body - divided into two parts, upper and lower bodys, the lower external body serves as a support to the work surface and the upper external body of the fume hood, built in steel tube and panels made of sheet steel coated with epoxy-polyester resins with high chemical resistance. The upper external body serves as support for the inner lining of the hood, made up of panels made of sheet steel coated with epoxy-polyester resins of high chemical resistance. In the back space of the bench there is a technical zone where connections to the water, sewage, gas and electricity networks are made;
- Inner lining - Made of various materials with high chemical resistance. It has aerodynamic elements to optimize the flow of air during the extraction.
- Working surface - Made of various materials with high chemical resistance.
- Window - Vertical opening (guillotine type), allows the user to handle dangerous substances inside the hood safely.

Our Standard models are supplied with a simple constant air volume (CAV) controller, and the services are a gas valve, and a water valve with the respective controls, it is also provided with a sink on the worktop, as well as, two or three IP55 single-phase sockets. The controls and electrical plugs are placed on the front panels located under the worktop of the fume hood, the nozzles being provided inside the hood on the worktop or panels on the sides of the inner lining. There is also the possibility to be supplied with the controls and sockets in the side towers.

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Air flow

Fluxo de ar
Flux d'air

The Standard Fume Hood comes with constant air volume (CAV), however can be equipped with VAV system (variable air volume), which improves the performance of the Fume Hood as well as its overall energy consumption.

In the same way the Fume Hood can also be equipped with automatic window closing system, or with the SCAT chemical waste collection system.

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

Standard services

Serviços standard
Services standard

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
One cold water valves and one natural gas valve, with the respective taps.

Uma válvulas de água fria e uma válvula de gás natural, com as respectivas torneiras.
Une vanne d'eau froide et une vanne de gaz naturel, avec les robinets correspondants.

One Drip Cup undermounted on the worktop.

Uma piletta fixa na superfície de trabalho.
Une bénitier fixé sur le plan de travail.

Electrical information*

Informação eléctrica*
Informations électriques*

The ENS Fume Hood is a class I equipment, so all accessible metal parts are connected to the earth terminal, which must be connected to the earth protective circuit of the building. The terminal strips for electrically connecting the Fume Hood are located in the electrical cabinet on the control panels.

The equipment has all necessary magnetothermic protections, namely:

- Thermal protection of the fan (s);
- Protection circuit breaker of the circuit of the sockets;
- Circuit breaker protection for lighting system;
- Protection circuit breaker for the control circuit.

The installation site must have a circuit-breaker to ensure that in the case of other equipment connected to the circuit-breaker, these shall not be without power in case of maintenance of the Fume Hood.

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


STEEL SHEET
CHAPA DE AÇO
TÔLE D'ACIER

STEEL SHEET
CHAPA DE AÇO
TÔLE D'ACIER

COMPACT LAMINATE
COMPACTO FENÓLICO
PHENOLIC COMPACT

1200x850x2450

103-0010

1500x850x2450

103-0011

1800x850x2450

103-0012

STEEL SHEET
CHAPA DE AÇO
TÔLE D'ACIER

STEEL SHEET
CHAPA DE AÇO
TÔLE D'ACIER

TECHNICAL CERAMIC
GRÉS CÉRÉMICO
CÉRAMIQUE TECHNIQUE

1200x850x2450

103-0004

1500x850x2450

103-0005

1800x850x2450

103-0006

STEEL SHEET
CHAPA DE AÇO
TÔLE D'ACIER

POLYPROPYLENE
POLIPROPILENO
POLYPROPYLENE

TECHNICAL CERAMIC
GRÉS CÉRÉMICO
CÉRAMIQUE TECHNIQUE

1200x850x2450

103-0034

1500x850x2450

103-0035

1800x850x2450

103-0036

STEEL SHEET
CHAPA DE AÇO
TÔLE D'ACIER

TECHNICAL CERAMIC
GRÉS CÉRÉMICO
CÉRAMIQUE TECHNIQUE

TECHNICAL CERAMIC
GRÉS CÉRÉMICO
CÉRAMIQUE TECHNIQUE

1200x850x2450

103-0021

1500x850x2450

103-0022

1800x850x2450

103-0023

L x D x H*

1200x850x2450

103-0031

1500x850x2450

103-0032

1800x850x2450

103-0033

STEEL SHEET
CHAPA DE AÇO
TÔLE D'ACIER

POLYPROPYLENE
POLIPROPILENO
POLYPROPYLENE

POLYPROPYLENE
POLIPROPILENO
POLYPROPYLENE

1200x850x2450

103-0031

1500x850x2450

103-0032

1800x850x2450

103-0033

L x D x H*

1200x850x2450

103-0037

1500x850x2450

103-0038

1800x850x2450

103-0039

STAINLESS STEEL
AÇO INOX
ACIER INoxyDABLE

STAINLESS STEEL
AÇO INOX
ACIER INoxyDABLE

STAINLESS STEEL
AÇO INOX
ACIER INoxyDABLE




*Length x Depth x Height
*Comprimento x Profundidade x Altura
*Longueur x Profondeur x Hauteur

STANDARD
STANDARD
STANDARD

STRONG ACIDS
ÁCIDOS FORTES
ACIDES FORTS

RADIOISOTOPE
RADIOSÓTOPO
RADIO-ISOTOPE

Fume Hood Indlab ENS Standard for a generalized use of chemical material, with the exterior in cold-rolled steel sheet and the inner lining and worktop in different materials. For a more intensive use of chemical materials it is recommended that the interior and the worktop be in ceramic stoneware, for a better performance of the fume hood over time.

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This ENS Fume Hood has been developed for acid chemical materials, with exterior in cold-rolled steel sheet, interior and work surface in stainless steel. It is a fume hood advised for those who intend to work with strong acids and very corrosive materials, ensuring stability and safety when used. However due to the nature of the polypropylene, it is not advisable to use material or equipment with high temperatures, since it will damage the polypropylene.

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This ENS Fume Hood was developed specifically for radioisotope chemical materials, with exterior, interior and work surface in stainless steel. It is a fume hood advised for those who work in high temperature environments and who need to be in constant cleaning, however, due to the nature of the material, it is not advised to work with some acids.

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