

EXTRACTION SYSTEMS

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Industrial laborum

CAV SYSTEM

The CAV system, (constant air volume), is the simplest system for extraction in fume hoods and other extraction systems, and as the name indicates is a system that maintains the volume of air extracted, and there may be variation in facial velocities.

For a system that uses a speed and volume of constant extraction, makes it necessary in the case of the fume hoods to resort to bypass so that there is no variation of the velocity to the face of the fume hood, less than 0,3m/s or greater than 0,5m/s. The Bypass thus allows the user to open and close the window without varying the speed and staying within the safety standards.

This system presents a reasonable level of safety but at economic levels becomes very expensive, since it makes extraction volumes that may be unnecessary if the window of the hood is closed, and a higher volume of extraction then means a volume of insufflation superior, all this afterwards is

translated to energy consumptions that in this type of system can be elevated.

As this extraction system presents constant air volumes, the HVAC system of the room where it is being inserted can be also simple, since to maintain the temperature and the atmospheric pressure within the parameters does not need to make great adjustments to what is normalized.

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
SYSTÈME CAV


Constant Air Volume
Volume de Ar Constante
Volume d'Air Constant

The CAV (Constant Air Volume) system is a constant air extraction system where the volume of air drawn by the fume hood is always constant regardless of the window opening. With the variation of the height of the window the speed of passage of air varies, for example, at a 500mm opening the facial speed is of 0.3 to 0.5 m / s, when closing the window the speed will increase, creating turbulences and noises. It is the most basic system, the least expensive in the acquisition and the most expensive on the usage.

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
Industrial Laborum features two types of controllers. The simple controller, which features only two buttons: one for turning on / off ventilation and the other for turning on / off the luminaire. The second CAV controller has a button to turn on / off the ventilation, another button to turn on / off the luminaire and has a screen that allows the visualization of the speed of air passage. This second presents an audible alarm for the air extraction failure. Both controllers will be implemented and secured outside the Fume Hood in one of the uprights or service panels below the work surface.

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Buttons
Botões
Boutons

ON/OFF Ventilation
Ventilação ON/OFF
Ventilation ON/OFF

ON/OFF Light
Luz ON/OFF
Lumière ON/OFF

Lights
Luces
Lumières

ON/OFF Ventilation
Ventilação ON/OFF
Ventilation ON/OFF

ON/OFF Light
Luz ON/OFF
Lumière ON/OFF

VAV SYSTEM

SISTEMA VAV
Système VAV

VAV SYSTEM

The VAV system (variable air volume) is the most advised system for an intensive and safe use of a Fume hood. With the volume variation there is a velocity control that the CAV system does not guarantees, this means that the user whatever the environment is always safe.

The system uses constant velocity readings, as soon as it is low, the sensor sends the information to the central control, which in turn passes the information to the actuator that adjusts the opening of the damper so that the speed returns to the desired levels.

In addition to the clear and obvious safety benefits of this system, the benefits are also economical. Since extractions and volumes of air are always adapted to the needs and there are no unnecessary consumptions, also the room HVAC must be a system with variable air volume, this way, we have an optimization of the consumptions, never consuming unnecessary.

Although it is a system that allows savings on a daily basis, it is a much more expensive system to buy than the CAV, however the savings in use ends up compensating the user a lot. In case of using a VAV system also the room HVAC will have to be thought of as VAV, so that the room parameters are met for your own safety.

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Variable Air Volume

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Volume d'Air Variable

As the name indicates, the VAV system it's a system that varies the volume of air extracted depending on the opening of the window, and thus we have a stable face velocity. This implies that the safety of the user is guaranteed with the window open or closed. The VAV system comes with a speed reader that makes a constant reading of the extraction values and communicates to the central system which in turn gives information to the actuator that opens or closes the damper, thus controlling the extraction volume and the facial velocity. This is the most advisable system for a laboratory, since it always guarantees the safety of the user and allows a better saving of energy costs.

As an extra safety measure when the window exceeds 500mm height or when conditions are not ideal for using the fume hood, the interior light changes from pure white to red, thereby alerting users, in a very visible way, that they should not be used under the conditions it is at that exact moment.

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



ILI VAV com variador de velocidade

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Cod. 120-0010

ILI VAV with monitoring & automatic sash

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ALARM
Red alarm light
Luz de alarme vermelha
Voyant d'alarme rouge

Controllers

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Industrial Laborum presents two types of controllers; the ILI controller which is the base VAV, the most economical version in the purchase, but nevertheless very different in the safety levels, compared to the premium systems. The second version of controller, are the premium versions, that features systems of different brands of renown in extraction systems in the world, they all have approximately the same characteristics, differing only in some aspects in the controller. The choice of VAV system should not be only monetary, as the system will be introduced into the internal HVAC system. We must take into account that the brand we are using in the HVAC, and try to make the brand of the VAV equal, since the compatibility systems will be much higher.

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SAFETY & ECONOMY

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


Industrial Laborum
lab integrated systems

Industrial Laborum

Industrial Laborum

Industrial Laborum

COMPARISON

COMPARAÇÃO

COMPARAISON

CAV

CAV

CAV

CAV

The CAV system is the cheapest at the time of purchase, but in the long time it becomes the most expensive. Since this does not adapt to the opening of the window the volume will have that air will have to be calculated for the maximum opening, this implies a high volume of extracted air. Thus the volume of air replacement in the room will be high, as well as, the temperature and pressure control in the room will have much higher costs compared to other systems.

The extraction made by the Fume Hoods with CAV system is an extraction that does not necessarily imply an HVAC of the room also governed by a CAV system; however, an HVAC with VAV system is most advised, since the control of the area where the inserts are as accurate as possible. This also enables the user to use only a hood without a sudden rise in temperature or an increase in pressure, which could put the user in danger.

All these variants later become high laboratory maintenance cost. It is thus summarized that the CAV system is not the most suitable for large laboratories, and that although it is cheaper in the short term, in the long term it becomes the most expensive.

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

VAV

VAV

Of all traditional extraction systems the most recommended is the VAV system. A system that makes constant measurements of the velocities inside the hood, thus being able to adapt the volume of extracted air. This adaptation of the system to the conditions allows for substantial energy savings and improved user safety, although this is a more expensive system than the CAV.

With a VAV extraction system the most advised is that the laboratory HVAC is also with VAV system, since thus the accuracy of the pressure, temperature and all environmental parameters of the room will be guaranteed and thus maintained.

The end result of this initial investment (higher priced systems) is that the user safety is improved and the cost of maintaining the hottes in use is much lower compared to the CAV system. This is the recommended system for large laboratories and for areas where atmospheric pressure and temperature control will need to be accurate.

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

GFH

GFH

The GFH system (Green Fumehood), is a filtration extraction system, developed from the latest technology in Erlab filters, the Neutrodine filters, a filter that allows the absorption of almost all chemical compounds and in greater quantity. This technology at economic levels implies a great investment in the beginning, but the system has an almost zero energy consumption, which in the end represents a very low maintenance cost.

In addition, the HVAC cost of the room is also reduced dramatically, since this system does not require extraction ducts and does not draw any air volume from the room, the HVAC system is reduced to temperature control and eventually air living room.

Summarizing the GFH system is a system that in the long term, even compared to the VAV system, ends up reducing the costs of use. However this system has to be well managed, and so when purchasing a Eco Smartflow Fume Hood will be made a report of maintenance costs of the hotte and compared to what it would cost if a traditional system is used, so that the acquisition of this equipment will ultimately compensate economically for the user.

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